



FINAL REPORT

Pretesting of molecular identification tests for *Dendrolimus sibiricus* (Chetverikov, 1908)

Contents

1. Introduction.....	1
2. Scope of pretesting	2
3. Test selection.....	2
4. Composition of the sample set	3
5. Specification of pretesting procedures	3
6. Results of pretesting.....	3
7. Database inventory for sequence records	4
8. Discussion.....	7
Appendix 1 - References.....	8
Appendix 2 – <i>In silico</i> data.....	9
Appendix 3 – Specifications and parameters for the molecular tests.....	14

1. Introduction

The European Reference Laboratory for Insects and Mites has to select, adapt or develop reliable identification tests for the phytosanitary insect and mite species that are relevant for the European Union (included in the Commission Delegated Regulation (EU) 2019/1702 and in the EURL for Insects and Mites working programmes). One of the tasks of the EURL is to validate available diagnostic protocols before recommending their use to the National Reference Laboratories of the European Union. Pretesting of available tests is necessary to select the most reliable ones for the validation study.

The Siberian conifer silk moth *Dendrolimus sibiricus* (Lepidoptera: Lasiocampidae) (Fig.1) is considered a fairly recently separated species from *Dendrolimus superans*, making their differentiation more difficult (this is also why the taxon is sometimes referred to as *Dendrolimus superans ssp. sibiricus*). Both are closely related to *D. pini*.



Fig. 1 *Dendrolimus sibiricus*, adult

D. sibiricus has thus far spread to China, Mongolia, North Korea and Russia, unlike *D. superans*, which is as of yet considered to be absent in Europe (only in north-east China, Japan and far eastern parts of Russia). So far, both species haven't been recorded in the European Union.

D. sibiricus is considered an European Union regulated species, listed among the EU quarantine pests (Annex II of the Commission Implementing Regulation (EU) 2019/2072) and among the EU priority pests (Commission Delegated Regulation (EU) 2019/1702).

More details are available in the EFSA Pest Survey Card on this pest (EFSA 2020).




Fig. 2 Current distribution of *D. sibiricus*

2. Scope of pretesting

The scope of this preliminary study was to identify molecular tests which are appropriate for the identification of *D. sibiricus*. Additionally, a database inventory for sequence records should shed light on the application possibilities of sequencing as identification methods.

3. Test selection

For this pest species no international standards are available yet, however they are currently under preparation.

The EPPO PM 7/129 (2) DNA barcoding as an identification tool for a number of regulated pests (EPPO, 2021) (Appendix 1), which includes tests for the DNA barcoding of arthropods in general was included in this pretesting.

According to Kononov *et al.* (2016), *D. sibiricus* can not be differentiated from *D. superans* and some *D. pini* genotypes sequencing the 5' end of the *COI* only.

ITS2 barcoding is mentioned in scientific literature as possible genetic locus for identification of *D. sibiricus* [Mikkola and Ståhls (2008), Kononov *et al.* (2016)].

Additionally, Mikkola and Ståhls (2008) give sequencing of the *COI* locus using CI-J-2183/ tI2-n-2014 (alias Jerry/Pat, called M5/3 by Kononov *et al.* 2016) primers as suitable for the identification of *D. sibiricus*.

For this pretesting study, ITS2 barcoding according to Kononov *et al.* (2016) was included, *COI* 3' end primers (CI-J-2183/ tI2-n-2014) were evaluated *in silico* only.

4. Composition of the sample set

For pretesting of tests indicated above, DNA was extracted from parts of single adult specimens (part of legs) of *Dendrolimus* spp. (Tab. 1).

Table 1: Details on single adult specimens for pretesting *Dendrolimus* spp. identification (sample set)

Sample	Species	Host plant	Origin	Source
2258/19	<i>D. sibiricus</i>	Pine	Russia	Kulinich, Lethmayer (AGES)
1590/20	<i>D. pini</i>	Pine	Estonia	Ermakovich (PHML)
1282/21	<i>D. sibiricus</i>	Pine	Russia (Omsk)	Filippov, Taddei (ANSES)
1283/21	<i>D. sibiricus</i>	Pine	Russia (Omsk)	Filippov, Taddei (ANSES)
1501/21	<i>D. sibiricus</i>	Pine	Russia (Omsk)	Filippov, Taddei (ANSES)
1502/21	<i>D. pini</i>	Pine	France	Vardon, Taddei (ANSES)
1503/21	<i>D. pini</i>	Pine	France	Vardon, Taddei (ANSES)
1504/21	<i>D. pini</i>	Pine	France	Vardon, Taddei (ANSES)

5. Specification of pretesting procedures

DNA extraction

For DNA extraction the DNeasy Blood and Tissue (Qiagen) on part of *Dendrolimus* spp. legs was used. Samples indicated in table 1 were tested undiluted and with at least one additional dilution (usually 1:20).

Molecular tests

Following tests for the identification of *D. sibiricus* were evaluated:

- *ITS2* sequencing according to Kononov *et al.* (2016)
- Included in EPPO PM 7/129 (2) DNA barcoding (EPPO, 2016): Barcoding for arthropods (Appendix 1): LepF/R primer set

In silico:

- 3' *COI* sequencing according to Mikkola and Ståhls (2008)

6. Results of pretesting

In silico data

Primer BLAST of sequencing primers included in the *in silico* pretesting of this study (LepF/R for the 5' *COI*, Pat/Jerry for the 3' *COI*, and ITS2A/B for the *ITS2* locus) indicated them as inclusive for a multitude of *Dendrolimus* species. Neighbor joining trees were constructed to evaluate the possibility of differentiating the closely related *D. superans*, *D. sibiricus* and *D. pini* species. For both the 3' and 5' ends of the *COI* locus, the separation seems was not possible. Only on the *ITS2* locus clustering on species level is visible (Appendix 2).

Testing of the sample set

With the sample set of this pretesting, the results of the *in silico* data as well as evaluations in literature could be confirmed. Whereas 5' *COI* barcoding was not able to distinguish between *D. sibiricus* and *D. superans*, correct identification of the *D. sibiricus* samples was possible on the *ITS2* locus (Tabl.2). Amplification of *D. pini* sequences on the *ITS2* locus proved difficult, however, but this species could in turn be correctly identified via 5' *COI* barcoding in most cases.

Table 2: Results of pretesting for *Dendrolimus* spp. according to assigned values

Assigned value of the samples	EPPO 2021, <i>COI</i> 5'	Kononov et al. 2016, <i>ITS2</i>
2258/19 <i>Dendrolimus sibiricus</i>	<i>D. superans</i>	<i>D. sibiricus</i>
1590/20 <i>Dendrolimus pini</i>	<i>D. superans</i>	<i>D. pini</i>
1282/21 <i>Dendrolimus sibiricus</i>	<i>D. superans</i>	<i>D. sibiricus</i>
1283/21 <i>Dendrolimus sibiricus</i>	<i>D. superans</i>	<i>D. sibiricus</i>
1501/21 <i>Dendrolimus sibiricus</i>	-	<i>D. sibiricus</i>
1502/21 <i>Dendrolimus pini</i>	<i>D. pini</i>	-
1503/21 <i>Dendrolimus pini</i>	<i>D. pini</i>	-
1504/21 <i>Dendrolimus pini</i>	<i>D. pini</i>	-

7. Database inventory for sequence records

DNA barcoding relies on PCR of predetermined marker genes (e.g. for the *COI* gene), DNA sequencing and comparison of those sequences to a database of reference sequences (Armstrong and Ball 2005). Applying barcoding for insect identification requires enough sequence records from the species within the genus for a reliable comparison. Not only the number, but also the genetic and geographic diversity of the records and the quality are potential issues that should be considered.

Three different databases (NCBI GenBank, Bold and EPPO Q-Bank) were consulted for the inventory. As search parameters the genus and species name and the gene locus (*COI* and synonyms) were used. In addition, the reliability of the records was checked and, if reasonable, the search was extended to other gene loci.

Results

In all three consulted databases sequence records for *D. sibiricus* on the *COI* locus are available, with both the 5' and 3' end being represented in GenBank and Bold. GenBank also offered sequences on the *ITS2* and other loci (e.g. *COII*, 5.8S rRNA) (query date October 2022).

Table 3: Number of sequence records per relevant gene for each database (data accessed October 2022).

Gene	GenBank	Q-Bank	Bold
<i>COI</i>	56	13	50
<i>COII</i>	16	0	0
<i>ITS2</i>	31	0	0

Detailed information

NCBI GenBank and Bold database hold sequence records on the *COI* locus for multiple species of the *Dendrolimus* genus (e.g. *D. punctatus*, *D. kikuchii*, *D. pini*, *D. sibiricus*, *D. superans*). In the EPPO Q-Bank 13 *COI* barcodes for *D. sibiricus* are deposited. 31 *ITS2* sequences are available in GenBank, in addition to several *COII* and 5.8S *rRNA* (included in the *ITS2* region) sequences.

Out of 184 *D. sibiricus* specimen found in the Bold database, 143 hold barcodes, 50 of which are public. The total *Dendrolimus* genus is credited with 1.047 specimens with barcodes out of 1.137 specimen records, belonging to 14 public species and an additional 18 non-public ones.

The geographic variation of the sequence records for *D. sibiricus* in Bold represents its known distribution, with the vast majority stemming from Russia and occasional specimens having been sampled in Mongolia and Kazakhstan. All records found in EPPO Q-Bank were larvae sampled in Russia.

Tree-based identification

To evaluate the species divergences within the genus, Neighbor Joining (NJ) trees of distance were constructed for the *ITS2* and *COI* loci using NCBI GenBank (max. seq. diff. of 0.75). Figure 3 confirmed the suitability of the *ITS2* locus for distinguishing multiple *Dendrolimus* species, whereas Figure 4 shows insufficient separation of *D. sibiricus* and *D. superans* on the *COI* locus.




Figure 3: Neighbor joining tree on the *ITS2* locus shows clustering of different *Dendrolimus* species in a way that allows clear separation.



Figure 4: Neighbor joining tree on the 5' COI locus shows insufficient separation of *D. sibiricus* and *D. superans*

Interestingly, many of the sequences for the 5' end of the COI locus deposited in the NCBI database and subsequently mined from there for the Bold database, included a gap of about 50bp in the middle of their 1371bp sequence (Fig. 5).

```

ACTCCTGGGTCTTTAATGGAGATGATCAAATTTATAATACTATTGTACACAGACATGCT
TTCAATCAATAATTTTTATAGTAATACCAATTAATAATGGAGGATTTGGTAATGGATTA
GTACCTTTAATATAGGGGCCCTGACATAGCATTCCACGAATAAATAATAAGATTT
TGATTAATACCCCTCTCTACCTTATTAATTTGAGGAAGAATTGTAGAAAGTGGAGCC
GGAACGTGGATGAACGTTTATCTCTCTATCTCTAATAATTTGCTCATAGAGGGAGATCA
GTTGATTTAGCTATTTTTCCCTTCATCTGCGGAATTTCACTAATTTAGGAAGCAAT
AATTTTATACAACAATTTAAACATAAAAATTAATAATATATCAATTTGATCAAACCT
TTATTTGTTGAGCTGATGGAATACAGCATTTTTATATTATATCACTACCAAGTCTT
GCCGGAGCAATTAATACTATTAACCTGATCGAAATTTAAATACATCAATCTTTGACCTT
GCTGGAGGAGGAGATCCCATTTTATATCAACATTTAATTTGATTTTTGGGNNNNNNNN
NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
AATAAAAAAGAAACCTTTGGTTGTTTAGGAATAAATTTAGCTATATTAGCAATTTGGATTA
TTAGGATTTATTTGATGAGCTCATCACATTTACAGTAGGATAGATATTGATACTCGA
GCTTATTTTACCTCTGCAACTATAAATTTGAGATACCAACAGGAATTTAAATTTTAGA
TGATTAGCAACCTCTCATGGAACACAAATAAATATAGCCCTTCCATCTTTGAAGATTA
GGATTTGATCTTATTTACTGTTGGAGGATTAACAGGATATCTTCCCAATTTCTCTT
ATCGATTTACCTCTCATGACTACTACTATGAGTAGCTCATTTTCAATGATGACTTTCA
ATAGGAGAGATTTGCTATTAAGGGGATTTATTCATTTGATACCCATTTTACTGGC
TTATCTATAAATCTTTTATATTAAAAATCAATTTTAAATATATTATTGGAGTAAAT
TTAACATTTTTCTCAACACTTTTAGGATAGCTGGAATACCTGACGATACTCTGAT
TACCCAGACTCATATATCTCATGAAATATTCTTCAATAGGTTCTTATATTCTTTA
TTAGGAGTTATAATAATTAATTTATTTGAGAATCAATAAATTAACCAAGCAATTAGA
ATCTTTACTTTAAACATAAGATCTTCTATTGAATGATACAAATTTACCT

```

Figure 5: Peculiarity recognised in several COI sequences of *D. sibiricus*.

8. Discussion

Molecular identification of *Dendrolimus sibiricus* (3' and 5' *COI* barcoding, *ITS2* sequencing) was pretested in this study. *ITS2* sequencing is crucial for the differentiation of *D. sibiricus* and *D. superans*, which is not possible on the 5' *COI* locus (typically used for barcoding). 3' *COI* sequencing was only tested *in silico* during this pretesting, but seems promising, providing enough resolution for these closely related species.

Additionally, the identification of the non-target *D. pini* is difficult using *ITS2* sequencing alone, as amplification of the sequence isn't always successful.

Based on this pretesting study, the use of more than one test seems necessary for accurate identification of *D. sibiricus*. Therefore, 5' *COI* barcoding combined with *ITS2* sequencing (Kononov *et al.* 2016) will be validated, in addition to 3' *COI* sequencing. The sample set for the validation will be extended to more specimens of both target and closely related non target Lepidoptera specimens.

Participating staff:

Experts/ Supervisors: Richard Gottsberger, Helga Reisenzein

Technical staff/ Operators: Claudia Heiss, Chiara Pohn

Date: February 8, 2023



Helga Reisenzein
EURL Deputy Director



Richard Gottsberger
Scientific Project Leader – Molecular Unit

Appendix 1 - References

EFSA (European Food Safety Authority), Wilstermann, A, Schrader, G, Kinkar, M, Vos, S, 2020. Pest survey card on *Dendrolimus sibiricus*. EFSA supporting publication 2020:EN-1779. 23 pp. doi: 10.2903/sp.efsa.2020.EN-1779.

EPPO (2021). EPPO standards PM 7/129 (2) DNA barcoding as an identification tool for a number of regulated pests. Bulletin OEPP/EPPO Bulletin, 51 (1): 100–143.

Hajibabaei, M., Janzen, D. H., Burns, J. M., Hallwachs, W., & Hebert, P. D. (2006). DNA barcodes distinguish species of tropical Lepidoptera. *Proceedings of the National Academy of Sciences*, 103(4), 968-971.

Kononov, A., Ustyantsev, K., Wang, B., Mastro, V. C., Fet, V., Blinov, A., & Baranchikov, Y. (2016). Genetic diversity among eight *Dendrolimus* species in Eurasia (Lepidoptera: Lasiocampidae) inferred from mitochondrial COI and COII, and nuclear ITS2 markers. *BMC genetics*, 17, 173-182.

Mikkola, K., & Ståhls, G. (2008). Morphological and molecular taxonomy of *Dendrolimus sibiricus* Chetverikov stat. rev. and allied lappet moths (Lepidoptera: Lasiocampidae), with description of a new species. *Entomologica Fennica*, 19(2), 65-85.

Simon, C., Frati, F., Beckenbach, A., Crespi, B., Liu, H. & Flook, P., 1994: Evolution, weighting and phylogenetic utility of mitochondrial gene-sequences and a compilation of conserved polymerase chain-reaction primers. — *Annals of the Entomological Society of America* 87: 651–701.

Appendix 2 – In silico data

Search set to *Dendrolimus* or *Dendrolimus sibiricus*, respectively.

5' COI sequencing:

EPPO (2021). EPPO standards PM 7/129 (2) DNA barcoding as an identification tool for a number of regulated pests. Bulletin OEPP/EPPO Bulletin, 51 (1): 100–143.

Hajibabaei, M., Janzen, D. H., Burns, J. M., Hallwachs, W., & Hebert, P. D. (2006). DNA barcodes distinguish species of tropical Lepidoptera. Proceedings of the National Academy of Sciences, 103(4), 968-971.

Fast Minimum Evolution tree for LepF:

Search set to *Dendrolimus*



3' COI sequencing:

Mikkola, K., & Ståhls, G. (2008). Morphological and molecular taxonomy of *Dendrolimus sibiricus* Chetverikov stat. rev. and allied lappet moths (Lepidoptera: Lasiocampidae), with description of a new species. *Entomologica Fennica*, 19(2), 65-85.

Simon, C., Frati, F., Beckenbach, A., Crespi, B., Liu, H. & Flook, P., 1994: Evolution, weighting and phylogenetic utility of mitochondrial gene-sequences and a compilation of conserved polymerase chain-reaction primers. — *Annals of the Entomological Society of America* 87: 651–701.

Fast Minimum Evolution trees for Primer forward: C1-J-2183 (alias Jerry, called M5 by Kononov *et al.* 2016)

5'-CAACATTTATTTTGGATTTTGG-3'

Search set to *Dendrolimus*

- moths | 17 leaves
- Dendrolimus superans isolate LBG100726.040 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate LBG100726.035 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate LBG100726.034 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate LBG100726.031 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate LBG100726.030 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate LBG100726.029 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate LBG100726.028 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate LBG100726.015 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate LBG100726.014 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus superans isolate JG110803.018 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus punctatus isolate DP cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii mitochondrion, complete genome
- Dendrolimus kikuchii isolate LY10 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate JI4 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate JI3 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate JL4 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate JL1 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate WYS11 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate WYS1 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate WYS14 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HY12 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate JI1 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate YS4 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate YS2 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HS22 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate QDH7 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HS15 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HS21 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate QDH4 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HS1 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate TP7 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HS6 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HS11 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HS22 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate HS22 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate QDH6 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate PE3 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate PE27 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus kikuchii isolate PE32 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus pini voucher 056/DENDSI/3B cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini voucher 056/DENDSI/3A cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus houii mitochondrion, complete genome
- Dendrolimus pini isolate FraXel cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate FraDig2 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate FraFDP21 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate BulgJ cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate BulgC cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate BulgM cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate TurkZolo cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate IotW11810 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate FraTS2 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate SwissMes3 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate Scot1 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate RusNovo18 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate SloveB cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate BulgB cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate BulgK cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate LithL2 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate SpaLLE cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate EngDD cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate GerWol31 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate RusRostov cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate FraCDM3 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate NorNW2 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate SweSE4 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate SpaInSpa10 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate RusZolo7 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate SwissWol Switzerland cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate GerKem2 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate SpaTer cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate NorNW3 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate LithZelva8 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate CroDub1 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate RusMosc2 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate DenBoto cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate LithZelva7 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate FraMLC1 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus pini isolate FraCP12 cytochrome c oxidase subunit I (COX1) gene, partial cds; mitochondrial
- Dendrolimus superans mitochondrion, complete genome
- icliQuery_21355

Search set to *D. sibiricus*



- Dendrolimus sibiricus isolate DS1 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus sibiricus isolate DS71 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus sibiricus isolate DS72 cytochrome oxidase subunit II (COII) gene, partial cds; mitochondrial
- Dendrolimus sibiricus isolate DS31 cytochrome oxidase subunit II (COII) gene, partial cds; mitochondrial
- Dendrolimus sibiricus voucher DS16188 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus sibiricus voucher DS14724 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus sibiricus voucher DS14706 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus sibiricus voucher DS14629 cytochrome oxidase subunit I (COI) gene, partial cds; mitochondrial
- Dendrolimus sibiricus mitochondrial partial coi gene for cytochrome c oxidase subunit I, specimen vouch...
- Dendrolimus sibiricus mitochondrial partial coi gene for cytochrome c oxidase subunit I, specimen vouch...
- Dendrolimus sibiricus voucher DS16188 5.8S ribosomal RNA gene, partial sequence; internal transcribe...
- Dendrolimus sibiricus voucher DS15586 5.8S ribosomal RNA gene, partial sequence; internal transcribe...
- Dendrolimus sibiricus voucher DS14708 5.8S ribosomal RNA gene, partial sequence; internal transcribe...
- Dendrolimus sibiricus voucher DS14631 5.8S ribosomal RNA gene, partial sequence; internal transcribe...
- Dendrolimus sibiricus voucher DS14629 5.8S ribosomal RNA gene, partial sequence; internal transcribe...
- Dendrolimus sibiricus isolate Novo11 internal transcribed spacer 2, partial sequence
- Dendrolimus sibiricus isolate DS32 internal transcribed spacer 2, partial sequence
- Dendrolimus sibiricus isolate DS31 internal transcribed spacer 2, partial sequence
- Dendrolimus sibiricus ITS2, specimen voucher MZH:KM5
- Dendrolimus sibiricus ITS2, specimen voucher MZH:KM31
- Dendrolimus sibiricus ITS2, specimen voucher MZH:KM1
- icliQuery_34641

Fast Minimum Evolution trees for Primer reverse: t12-n-3014 (alias Pat, called M3 by Kononov *et al.* 2016)

5'-TCCAATGCACTAATCTGCCATATTA-3'

Search set to *Dendrolimus*:

Search set to *Dendrolimus sibiricus*




ITS2 sequencing:

Kononov, A., Ustyantsev, K., Wang, B., Mastro, V. C., Fet, V., Blinov, A., & Baranchikov, Y. (2016). Genetic diversity among eight *Dendrolimus* species in Eurasia (Lepidoptera: Lasiocampidae) inferred from mitochondrial COI and COII, and nuclear ITS2 markers. *BMC genetics*, 17, 173-182.


Mikkola, K., & Ståhls, G. (2008). Morphological and molecular taxonomy of *Dendrolimus sibiricus* Chetverikov stat. rev. and allied lappet moths (Lepidoptera: Lasiocampidae), with description of a new species. *Entomologica Fennica*, 19(2), 65-85.

Fast Minimum Evolution trees for ITS2A:


Search set to *Dendrolimus*




Search set to *Dendrolimus sibiricus*



Fast Minimum Evolution tree for ITS2B:
Search set to *Dendrolimus*



Search set to *Dendrolimus sibiricus*



Appendix 3 – Specifications and parameters for the molecular tests

Specification of the PCR Assay 1 (COI Barcoding according to EPPO PM7/129(2))

Name of the primer incl. sequence, literature reference, fragment length in bp:

LepF: 5'- ATTCAACCAATCATAAAGATATTGG-3'

LepR: 5'- TAAACTTCTGGATGTCCAAAAAATCA-3'

Literature: Hajibabaei, M., Janzen, D. H., Burns, J. M., Hallwachs, W., & Hebert, P. D. (2006). DNA barcodes distinguish species of tropical Lepidoptera. *Proceedings of the National Academy of Sciences*, 103(4), 968-971.

Fragment length: 709bp

PCR - Parameters:

Thermocycler used: Biometra T3000 Thermal cycler

Mastermix: 5x HOT FIREPol® Master Mix, Solis Biodyne:

Composition:		Final concentration:
	Volume per reaction µl	
Water	6	
Mastermix	2	1x
Primer1:	0,5	0,5µM
Primer2:	0,5	0,5µM
∑	9	
DNA	1	

PCR conditions:

	°C	Duration (min., sec.)	Nr. of Cycles
Start	95	15 min	1
Denaturation	95	45 sec	5
Annealing	44	45 sec	
Extension	72	45 sec	
Denaturation	95	45 sec	35
Annealing	49	45 sec	
Extension	72	45 sec	
Final extension	72	7 min	1
Cooling	15	∞	

Specification of the PCR Assay 2 (ITS sequencing according to Kononov *et al.* (2016))

Name of the primer incl. sequence, literature reference, fragment length in bp:

ITS2A: 5'- TGTGAACTGCAGGACACAT-3'

ITS2B: 5'- TATGCTTAAATTGAGGGGGT-3'

Kononov, A., Ustyantsev, K., Wang, B., Mastro, V. C., Fet, V., Blinov, A., & Baranchikov, Y. (2016). Genetic diversity among eight *Dendrolimus* species in Eurasia (Lepidoptera: Lasiocampidae) inferred from mitochondrial COI and COII, and nuclear ITS2 markers. *BMC genetics*, 17, 173-182.

Fragment length: 504bp

PCR - Parameters:

Thermocycler used: Biometra T3000 Thermal cycler

Mastermix: 5x HOT FIREPol® Master Mix, Solis Biodyne:

Composition:		Final concentration:
	Volume per reaction μ l	
Water	6	
Mastermix	2	1x
Primer1:	0,5	0,5 μ M
Primer2:	0,5	0,5 μ M
Σ	9	
DNA	1	

PCR conditions:

	$^{\circ}$ C	Duration (min., sec.)	Nr. of Cycles
Start	95	15 min	1
Denaturation	95	30 sec	40
Annealing	52	30 sec	
Extension	72	30 sec	
Final extension	72	7 min	1
Cooling	15	∞	