



# Sanger Sequencing of the full-length rabies virus nucleoprotein gene

**Version 02**

September 2026



**ANSES - Nancy Laboratory for Rabies and Wildlife**  
**European Union Reference Laboratory for Rabies**

# CONTENTS

|                                                      |                      |
|------------------------------------------------------|----------------------|
| Contents.....                                        | 2                    |
| Foreword.....                                        | <a href="#">33</a>   |
| Introduction.....                                    | <a href="#">33</a>   |
| 1. Purpose and scope .....                           | <a href="#">44</a>   |
| 2. Reference documents .....                         | <a href="#">44</a>   |
| 3. Abbreviations .....                               | <a href="#">55</a>   |
| 4. Materials.....                                    | <a href="#">55</a>   |
| 4.1. Biological material used .....                  | <a href="#">55</a>   |
| 4.2. Reagents/material used .....                    | <a href="#">55</a>   |
| 4.3. Equipment .....                                 | <a href="#">66</a>   |
| 5. Operating Procedure.....                          | <a href="#">66</a>   |
| 5.1. OneStep RT-PCR.....                             | <a href="#">66</a>   |
| 5.1.1. Primers.....                                  | <a href="#">66</a>   |
| 5.1.2. Preparation of the master mix.....            | <a href="#">77</a>   |
| 5.1.3. Addition of template.....                     | <a href="#">77</a>   |
| 5.1.4. Setting up of thermocycler.....               | <a href="#">77</a>   |
| 5.2. Sequencing PCR.....                             | <a href="#">88</a>   |
| 5.2.1. Primers.....                                  | <a href="#">88</a>   |
| 5.2.2. Preparation of the master mix.....            | <a href="#">88</a>   |
| 5.2.3. Addition of template.....                     | <a href="#">99</a>   |
| 5.2.4. Setting up of conventional thermocycler ..... | <a href="#">99</a>   |
| 5.3. Electrophoresis .....                           | <a href="#">99</a>   |
| 5.4. Analysis of PCR products .....                  | <a href="#">1010</a> |
| 6. Results .....                                     | <a href="#">1010</a> |
| 6.1. Test validation .....                           | <a href="#">1010</a> |
| 6.2. Sequencing of PCR products .....                | <a href="#">1010</a> |
| 6.3. Analysis of DNA sequencing results .....        | <a href="#">1111</a> |

## **FOREWORD**

This protocol has been adapted by:

ANSES - Nancy Laboratory for Rabies and Wildlife

National Reference Laboratory and European Union Reference Laboratory for Rabies

Address: Technopole agricole et Vétérinaire, 54220 MALZEVILLE, France.

Version 2 : Update of reagent volumes included in the reaction mixtures

## **INTRODUCTION**

This Standard Operating Procedure (SOP) describes the Sanger sequencing of the full RABV nucleoprotein gene. Sanger sequencing also known as the "chain termination method" is based on the detection of labelled chain-terminating nucleotides that are incorporated by a DNA polymerase during the replication of the template submitted to the analysis. This method developed in 1977 by Frederick Sanger and his colleagues is the "first-generation" DNA sequencing method and has been extensively used for more than 40 years across the world. Despite the advantages of the next-generation sequencing (NGS), Sanger sequencing is still widely used, particularly when partial genome characterization is sufficient for virus typing or for sequencing of single genes and in some cases for the analysis of longer fragments (~1000 bp in length).

This protocol should not be considered transferable as it stands. It shall be adapted by each user for its working conditions and shall be validated with reference materials or samples of known status (as performed in the frame of rabies diagnosis proficiency testing by example).

All commercial references described in this SOP are given for information. Other equivalent reagents, or equipment could be used as far as it does not affect the results.

## **1. PURPOSE AND SCOPE**

This SOP describes the preparation of PCR products from positive samples for rabies virus (RABV) in order to characterise RABV variants. Many national and international service providers with competitive costs carry out routinely Sanger sequencing and usually achieved the DNA sequencing results at D+1. This SOP will therefore be limited to the preparation of PCR products for sending and sequencing by a service provider for the virus characterisation.

The SOP is based on a conventional (gel-based) PCR assay and the use of specific and sequencing primers. The primers used in RT-PCR are specific of RABV; for sequencing PCR, the primers used are modified RABV primers, combined with primers (M13uni-21; M13rev-29) conventionally used for Sanger sequencing. Sanger sequencing is carried out using the primers used in the sequencing PCR and internal primers specific to RABV.

This method is based on viral RNA extracted from organs and tissues.

The protocol described here is intended using the MasterCycler<sup>R</sup> thermal cycler (Eppendorf, France).

## **2. REFERENCE DOCUMENTS**

1. OIE. Chapter 3.1.18. Rabies (Infection with rabies virus and other Lyssaviruses). Manual of Diagnostic Tests and Vaccines for Terrestrial Animals 2023 Paris, France; World Organisation for Animal Health.
2. Darryl L. Davis, Edward P. O'Brien, and Catherine M. Bentzley. Analysis of the degradation of oligonucleotide strands during the freezing/thawing processes using MALDI-MS. Anal Chem 2000 Oct 15;72(20):5092-6.

### **3. ABBREVIATIONS**

|        |                                                 |
|--------|-------------------------------------------------|
| bp     | base paires                                     |
| RT-PCR | reverse transcriptase polymerase chain reaction |
| ICTV   | International Committee on Taxonomy of Viruses  |
| N      | Nucleoprotein gene                              |
| PCR    | PCR                                             |
| RABV   | Classical rabies virus                          |
| RNA    | Ribonucleic Acid                                |

### **4. MATERIALS**

#### **4.1. BIOLOGICAL MATERIAL USED**

- Viral RNA extracted from the supernatant of a clarified 10% brain tissue suspension or from brain samples stored in RNA Later
- Target positive controls:
  - Process sentinel (Positive control of process): Viral RNA extracted from a brain tissue suspension artificially contaminated by RABV at a concentration of 1-100x the LDMethod
  - Amplification sentinel (Positive control of RT-PCR): Synthetic RABV RNA (*in vitro* transcribed RNA) used at a concentration of 1-10x the LDPCR
- Non-target positive control:
  - Viral RNA extracted from an uninfected clarified brain suspension (i.e. mouse, birds, ...)
- Target negative control:
  - Viral RNA extracted from an uninfected clarified brain suspension (i.e. mouse, birds, ...)

#### **4.2. REAGENTS/MATERIAL USED<sup>1</sup>**

- RT-PCR kit: QIAGEN OneStep RT-PCR kit (Cat: 210212, Qiagen, France)
- Nested PCR: Platinum Taq DNA polymerase and its associated buffer (Cat: 10966034, Invitrogen, France)
- RNase Out 40 U/μl (Cat: 10777019, Invitrogen, France)
- dNTP 10mM (RT-PCR) (sequencing PCR)
- Ultra pure DNase-RNase free distilled water
- Electrophoresis: Loading dye buffer, agarose, nucleic acid gel staining

---

<sup>1</sup> All commercial references described in this document are given for information. Other equivalent reagents or equipment can be used. Each user must validate its protocol in its working conditions (reagents, equipment).

### 4.3. **EQUIPMENT**

- Dedicated and separated laboratory rooms or workstations for the preparation of primers/master mix, extraction of RNA, PCR amplification and post PCR
- PCR workstations with specific dedicated micropipettes
- Spin micro-centrifuges, refrigerated centrifuges
- Vortex
- Electrophoresis apparatus
- Conventional thermocycler

## 5. **OPERATING PROCEDURE**

The method described here is the protocol routinely used in the laboratory to characterise the variants of RABV using amplicons of 1,353 bp of length.

### 5.1. **ONESTEP RT-PCR**

The protocol described here is intended for use with the Eppendorf MasterCycler<sup>R</sup>.

#### 5.1.1. **PRIMERS**

The RT-PCR is performed with HPSF primers.

The molecular models used to **characterise variants of the rabies virus** are described below:

#### **Molecular models used:**

| Primer   | Sense | Sequence (5' 3')          | Position  | Function | Target     |
|----------|-------|---------------------------|-----------|----------|------------|
| JW12     | F     | ATG TAA CAC CYC TAC AAT G | 55-73     | RT-PCR   | Lyssavirus |
| PVN8-Gt1 | R     | AGT YTC TTC RGC CAT CTC   | 1585-1568 | RT-PCR   | RABV       |

*F: forward; R: reverse ; M13<sup>(1)</sup>: M13uni-21; M13<sup>(2)</sup>: M13rev-29*

All the masters stocks are at 100µM, aliquoted and stored at < -16°C. Working primers are at 20µM, aliquoted and stored at < -16°C.

Repeated freeze/thaw cycles should be avoided. In the EURL, the number of freezing / thawing is limited to 5 times in order to preserve the characteristics of primers.

### **5.1.2. PREPARATION OF THE MASTER MIX<sup>2</sup>**

1. Thaw all frozen reagents (i.e. Qiagen One step RT-PCR, primers, molecular grade water, dNTPs). Mix the individual solutions and place them on ice until the addition of the enzyme mix.
2. Prepare the master mix as follows:

| Master mix                             | Final concentration | Volume (in $\mu\text{L}$ ) for 1 tube |
|----------------------------------------|---------------------|---------------------------------------|
| DNase/RNase free water                 |                     | 4.3                                   |
| QIAGEN OneStep RT-PCR (5X)             | 1X                  | 3                                     |
| dNTPs (10mM)                           | 0.4mM               | 0.6                                   |
| JW12 (20 $\mu\text{M}$ )               | 0.7 $\mu\text{M}$   | 0.5                                   |
| PVN8-Gt1 (20 $\mu\text{M}$ )           | 0.7 $\mu\text{M}$   | 0.5                                   |
| Qiagen RT-PCR Ez mix                   |                     | 0.6                                   |
| RNase Inhibitor (40 U/ $\mu\text{l}$ ) | 20 U                | 0.5                                   |

3. Mix the reaction thoroughly.
4. Dispense a volume of 10 $\mu\text{L}$  of master mix into PCR tubes.

### **5.1.3. ADDITION OF TEMPLATE**

5. Add 5 $\mu\text{L}$  of RNA samples/control RNA(s) to the individual PCR tubes containing the master mix and close the tubes.
6. Transfer the PCR tubes to the thermocycler for thermal cycling.

**For each run, include positive and negative controls** as follows:

- A **PCR negative control** = ultrapure water instead of sample to validate the no-contamination of the PCR.
- A **PCR positive control** = known target RNA previously extracted and included to ensure the good running of PCR.

### **5.1.4. SETTING UP OF THERMOCYCLER<sup>3</sup>**

7. Load the samples into the machine.
8. Program the thermocycler as follows:

---

<sup>2</sup> There are numerous available commercially reagents dedicated for RT-PCR and PCR. Optimisation of PCR (master mix formulation and cycling parameters) should be carried out by each user according to the enzyme mix used and the thermocycler used under his own working conditions.

| STEP                  | Temperature | Time   | Number of cycles |
|-----------------------|-------------|--------|------------------|
| Reverse transcription | 50°C        | 30 min | 1                |
| Polymerase activation | 95°C        | 15 min | 1                |
| Denaturation          | 94°C        | 30 sec | 16x              |
| Annealing             | 57°C- 51°C* | 30 sec |                  |
| Extension             | 72°C        | 2 min  |                  |
| Denaturation          | 94°C        | 30 sec | 29x              |
| Annealing             | 55°C        | 30 sec |                  |
| Extension             | 72°C        | 2 min  |                  |
| Final extension       | 72°C        | 10 min | 1                |

\*Touch-down PCR : decrease of 0.4°C /cycle

**Storage:** RT-PCR products are used to implement the sequencing PCR. They are stored at 5±3°C if used during the day or at <-16±3°C for longer storage.

## 5.2. SEQUENCING PCR

The protocol described here is intended for use with the Eppendorf MasterCycler<sup>R</sup>.

### 5.2.1. PRIMERS

**Molecular models used:**

| Primer                            | Sense | Sequence (5' 3')                                     | Position | Function   | Target |
|-----------------------------------|-------|------------------------------------------------------|----------|------------|--------|
| M13 <sup>(1)</sup> uni-21_JW12    | F     | TGT AAA ACG ACG GCC AGT ATG<br>TAA CAC CYC TAC AAT G | /        | Second PCR | RABV   |
| M13 <sup>(2)</sup> rev-29_PVN8bis | R     | CAG GAA ACA GCT ATG ACC TTG<br>CTC ATA YTT GGG       | /        | Second PCR | RABV   |

All the masters stocks (dedicated to the sequencing RT-PCR) are at 100µM, aliquoted and stored at < -16°C. Working primers are at 20µM, aliquoted and stored at < -16°C.

Repeated freeze/thaw cycles should be avoided. In the EURL, the number of freezing / thawing is limited to 5 times in order to preserve the characteristics of primers.

### 5.2.2. PREPARATION OF THE MASTER MIX

9. Mix the individual solutions and place them on ice.
10. Prepare the master mix as follows:

<sup>3</sup> The cycling parameters have been validated on the conventional thermocycler Eppendorf Mastercycler Gradient. Thermal profiles may be subject to optimization depending on used PCR equipment and kits.

| Master mix                               | Final concentration | Volume (in $\mu\text{L}$ ) for 1 tube |
|------------------------------------------|---------------------|---------------------------------------|
| DNase/RNase free water                   |                     | 14.05                                 |
| PCR Buffer without $\text{MgCl}_2$ (10X) | 1X                  | 2                                     |
| $\text{MgCl}_2$ 50mM                     | 2.5mM               | 1                                     |
| dNTPs 10mM                               | 0.1mM               | 0.20                                  |
| M13-JW12 (20 $\mu\text{M}$ )             | 0.25 $\mu\text{M}$  | 0.25                                  |
| M13-PVN8bis (20 $\mu\text{M}$ )          | 0.25 $\mu\text{M}$  | 0.25                                  |
| Taq Platinum (5U/ $\mu\text{l}$ )        | 1.25U               | 0.25                                  |

11. Mix the reaction thoroughly.
12. Dispense a volume of 18 $\mu\text{L}$  of master mix into PCR tubes.

### 5.2.3. ADDITION OF TEMPLATE

13. Add 2 $\mu\text{L}$  of the RT-PCR products to the individual PCR tubes containing the master mix and close the tubes.
14. Transfer the PCR tubes to the thermocycler for thermal cycling.

### 5.2.4. SETTING UP OF THERMOCYCLER

15. Load the samples into the machine.
16. Program the thermocycler as follows:

| STEP                         | Temperature | Time   | Number of cycles |
|------------------------------|-------------|--------|------------------|
| <b>Polymerase activation</b> | 94°C        | 2 min  | 1                |
| <b>Denaturation</b>          | 94°C        | 30 sec | 45x              |
| <b>Annealing</b>             | 48°C        | 30 sec |                  |
| <b>Extension</b>             | 72°C        | 2 min  |                  |
| <b>Final extension</b>       | 72°C        | 7 min  | 1                |

17. Transfer the PCR tubes to a room dedicated to electrophoresis.

**Storage:** PCR products are stored at a temperature of  $5\pm 3^\circ\text{C}$  until the realization of the electrophoresis then at  $-16^\circ\text{C} \pm 3^\circ\text{C}$  before sending to the sequencing provider.

## 5.3. ELECTROPHORESIS

5  $\mu\text{L}$  of PCR products are added to 1  $\mu\text{L}$  of loading dye buffer (loading buffer is 6x concentrated). The mixture is then applied to a 1.5% agarose gel, to which an intercalant<sup>4</sup> has been added.

<sup>4</sup> Ethidium bromide has been used for decades to stain nucleic acids in agarose gels under UV light. However, because ethidium bromide has been shown to be mutagenic and potentially carcinogen, the safer SYBR-based options have been preferred by many users over the use of ethidium bromide.

The nucleic acid gel stain used in the EURL is the SYBR safe at a final concentration of 1:10.000 (i.e. 1µl SYBR Safe per 10ml of agarose gel).

A well is reserved for the deposition of a DNA size marker in order to estimate the size of amplicons.

#### **5.4. ANALYSIS OF PCR PRODUCTS**

**The result is positive when a DNA fragment is detected at the expected size of 1520 bp (M13-JW12/M13-PVN8bis).**

**The result is negative when no DNA fragments are visible for the nested PCR.**

### **6. RESULTS**

#### **6.1. TEST VALIDATION**

The tests are validated as follows:

- if the positive control = positive.
- if the negative control = negative.

#### **6.2. SEQUENCING OF PCR PRODUCTS**

Once the test is validated, the PCR products are sent to a service provider for Sanger sequencing following the recommendations of the service provider.

Sequencing should be carried out in each direction to minimise the risk of errors in the sequences. Use the sequencing primers (forward M13uni-21\_JW12 and reverse M13rev-29\_PVN8bis) as well as the four internal primers described below. The PCR products should be purified prior to sequencing. Usually, the service provider carry out purification of PCR products (usually, free of charge).

## Molecular models used:

| Primer      | Sense | Sequence (5' 3')              | Function   | Target |
|-------------|-------|-------------------------------|------------|--------|
| M13uni-21   | F     | TGT AAA ACG ACG GCC AGT       | Sequencing | RABV   |
| M13rev-29   | R     | CAG GAA ACA GCT ATG ACC       |            |        |
| JW12_N1500  | F     | ATG TAA CAC CYC TAC AAT G     |            |        |
| JW662-684_F | F     | CAC GGT TGT TAC TGC TTA TG    |            |        |
| JW938-959_F | F     | GTG TTC AAT CTC ATT CAC TTT G |            |        |
| JW938-959_R | R     | CAA AGT GAA TGA GAT TGA ACA C |            |        |

F: forward; R: reverse ; M13<sup>(1)</sup>: M13uni-21; M13<sup>(2)</sup>: M13rev-29

DNA concentration and volume should respect the service provider's recommendations. DNA quantification can be estimated via agarose gel or a photometer to ensure accurate results.

An example of volume and sample concentration recommended by a service provider is shown for information below:

- Plasmid DNA with a product length < 30 kbp:
  - o Recommended sample concentration: 50-100 ng/μL,
  - o Volume of 15 μL for 1-4 reactions
- Purified PCR products:
  - o 150-300bp: 1 ng/μL
  - o 300-1000bp: 5 ng/μL
  - o 1000-3000bp: 10 ng/μL
- Unpurified products:
  - o 150-300bp: 4 ng/μL
  - o 300-1000bp: 10 ng/μL
  - o 1000-3000bp: 20 ng/μL

### 6.3. ANALYSIS OF DNA SEQUENCING RESULTS

The raw data can be provided in different formats depending on the model of sequencer used. Usually, the service provider provides a chromatogram, a sequence in the format ".fasta" and a sequence in the format ".ab1" or ".seq" from which we can analyse and derive the sequence.

In the EURL, the analysis of data and the generation of consensus sequence is performed using the Geneious Prime software.

Numerous software packages exist for the analysis of sequences as well as the generation of consensus sequences (e.g. Vector NTI from Invitrogen, Seqman in DNASTAR Lasergene).