

# STANDARD OPERATING PROCEDURE FMD VP1 AMPLIFICATION

## 1. PURPOSE AND SCOPE OF APPLICATION

This test is used for the virological diagnosis of foot-and-mouth disease (FMD) in cases of clinical suspicion (as a second-line test) or for epidemiological surveillance in the context of a diagnostic or research project. This test allows the amplification of the VP1 coding region of the FMD virus genome. The sequence analysis of this region allows to perform phylogenetic study and the characterization of the serotype/topotype/lineage of the FMD virus.

## 2. BIBLIOGRAPHIC REFERENCES

Analytical method recommended for molecular epidemiological study by the World Organisation for the Animals Health (WOAH) in the Manual of Diagnostic Tests and Vaccines for Terrestrial Animals (English version May 2022).

### *Bibliographic references:*

- 1) Ayelet G, Mahapatra M, Gelaye E, Egziabher BG, Rufeal T, Sahle M, Ferris NP, Wadsworth J, Hutchings GH, and Knowles NJ. Genetic Characterization of Foot-and-Mouth Disease Viruses, Ethiopia, 1981-2007. Technical appendix (Table 2) Emerging Infectious Diseases • www.cdc.gov/eid • Vol. 15, No. 9, September 2009. doi: 10.3201/eid1509.090091
- 2) Knowles NJ, Wadsworth J, Bachanek-Bankowska K, King DP. VP1 sequencing protocol for foot and mouth disease virus molecular epidemiology. Rev Sci Tech. 2016 Dec; 35 (3):741-755. doi: 10.20506/rst.35.3.2565
- 3) Valarcher JF, Knowles NJ, Zakharov V, Scherbakov A, Zhang Z, Shang YJ, Liu ZX, Liu XT, Sanyal A, Hemadri D, Tosh C, Rasool TJ, Pattnaik B, Schumann KR, Beckham TR, Linchongsubongkoch W, Ferris NP, Roeder PL, Paton DJ. Multiple origins of foot-and-mouth disease virus serotype Asia 1 outbreaks, 2003-2007. Emerg Infect Dis. 2009 Jul; 15 (7):1046-51. Appendix.
- 4) EL-Bayoumy MK, Abdelrahman KA, Allam AM, Farag TK, Abou-Zeina Hala AA and Kutkat MA. Molecular Characterization of Foot-and-Mouth Disease Virus collected from Al-Fayoum and Beni- Suef Governorates in Egypt. Global Veterinaria 13 (5): 828-835, 2014. DOI: 10.5829/idosi.gv.2014.13.05.86188.
- 5) Wadsworth J, Knowles NJ, Swabey KG, Stirling JM, Statham RJ, Li Y, Hutchings GH, Ferris NP, Paton DJ. Recent spread of new strains of foot-and-mouth disease virus type A in the Middle East and North Africa. Appendix 15 ([http://www.fao.org/ag/againfo/commissions/docs/research\\_group/paphos/app15.pdf](http://www.fao.org/ag/againfo/commissions/docs/research_group/paphos/app15.pdf))
- 6) Knowles NJ and Samuel AR. RT-PCR and Sequencing Protocols for the Molecular Epidemiology of Exotic Virus Diseases of Animals. 1998. ([http://www.wrlfmd.org/fmd\\_genotyping/genotyping\\_protocol.pdf](http://www.wrlfmd.org/fmd_genotyping/genotyping_protocol.pdf)).

## 3. PRINCIPLE OF THE METHOD

After RNA extraction, the test is first based on reverse transcription (RT) of the FMD virus RNA into complementary DNA (cDNA). The cDNA obtained is then amplified by PCR using a DNA polymerase that

uses primers targeting the VP1 region that varies according to the virus serotype, thus allowing the specific amplification of the 7 serotypes of FMDV. The two enzymatic reactions are performed successively in the same tube (one-step RT-PCR).

Due to the high rate of mutations during the replication of RNA viruses, it is recommended to use simultaneously at least 2 of the 3 pairs of primers specific to the serotype searched

## 4. REAGENTS AND PRIMERS

### 4.1 Reagents

One-step RT-PCR kit, QIAGEN, ref 210212, stored at – 21 °C (± 5 °C).

Ultra pure agarose, Invitrogen (or equivalent)

Loading buffer 6x, Biolabs ref B7021S (or equivalent)

TAE 1X prepared from concentrated TAE 50X, Euromedex ref EU0201-B (or equivalent)

Molecular weight marker, Marker VI, Roche, ref 11062590001 (or equivalent)

Ethidium bromide, BIOPROBE, ref. ETBC01, diluted at 1/10 000 in distilled water (or equivalent).

### 4.2 Primers

All primer stocks are at 100 µM, aliquoted by 100 µL / tube maximum and stored at -21°C (± 5°C). As soon as the penultimate aliquot is used up, a new primer tube should be ordered (alarm stock = 2 aliquots).

The primers used for the amplification of the genome fragment to be sequenced are described in tables 1 and 2.

| Serotype      | Sens | Name       | Sequence (5'-3')           | Region | Size expected | Alignment  | Ref.    |
|---------------|------|------------|----------------------------|--------|---------------|------------|---------|
| O             | +    | O-1C-244F  | GCAGCAAACACATGTCAAACACCTT  | VP3    | 1180 bp       | AY593823.1 | 1, 2    |
| O             | +    | O-1C-272F  | TBGCRGGNCTYGCCAGTACTAC     | VP3    | 1152 bp       | AY593823.1 | 1, 2    |
| O             | +    | O-1C-283F  | GCCCAGTACTACACACAGTACAG    | VP3    | 1141 bp       | AY593823.1 | 1, 2    |
| A             | +    | A-1C-562F  | TACCAAATTACACACGGGAA       | VP3    | 874 bp        | AY593762.1 | 1, 2    |
| A             | +    | A-1C-612F  | TAGCGCCGGCAAAGACTTTGA      | VP3    | 824 bp        | AY593762.1 | 1, 2    |
| C             | +    | C-1C-536F  | TACAGGGATGGGTCTGTGTGTACC   | VP3    | 945 bp        | AY593804.1 | 1, 2    |
| C             | +    | C-1C-616F  | AAAGACTTTGAGCTCCGGCTACC    | VP3    | 865 bp        | AY593804.1 | 1, 2    |
| Asia1         | +    | As1-1C505F | TACACTGCTTCTGACGTGGC       | VP3    | 888 bp        | JF739177.1 | 3       |
| Asia1         | +    | As1-1C530F | CCACRAGTGTGCARGGATGGGT     | VP3    | 863 bp        | JF739177.1 | 3       |
| O, A, C, Asia | -    | EUR-2B52R  | GACATGTCCTCCTGCATCTGGTTGAT | 2B     | /             | /          | 1, 2, 3 |

**Table1:** Recommended primers for amplification and sequencing of the VP1 region of serotypes O, A, C and Asia.

| Serotype | Sens | Name          | Sequence (5'-3')         | Region | Size expected | Alignment | Ref. |
|----------|------|---------------|--------------------------|--------|---------------|-----------|------|
| SAT1     | +    | SAT1-1C-559F  | GTGTATCAGATCACAGACACACA  | VP3    | 1042bp        | AY593845  | 1, 2 |
| SAT1     | +    | SAT1-1U-OSF   | GTGTACCAGATCACTGACAC     | VP3    | 1042bp        | AY593845  | 1, 2 |
| SAT1     | +    | SAT1-P1-1228F | AACCTGCACTTCATGTACAC     | VP3    | 1285bp        | AY593845  | 2    |
| SAT2     | +    | SAT2-1C-445F  | TGGGACACMGGIYTGAACTC     | VP3    | 1125bp        | AF540910  | 1, 2 |
| SAT2     | +    | SAT2-P1-1223F | TGAACTACCACTTCATGTACACAG | VP3    | 1259bp        | AF540910  | 1, 2 |
| SAT2     | +    | SAT2-VP3-AB F | CACTGCTACCACTCRGAGTG     | VP3    | 1143bp        | AF540910  | 1    |

|          |   |               |                       |     |        |          |      |
|----------|---|---------------|-----------------------|-----|--------|----------|------|
| SAT3     | + | SAT3-1C-1222F | AATCTGCATTTTCATGTACAC | VP3 | 1128bp | AY593850 | 2    |
| SAT3     | + | SAT3-1C-559F  | CTGTACCAAATYACAGACAC  | VP3 | 1280bp | AY593850 | 2    |
| SAT1-2-3 | - | SAT2B208R     | ACAGCGGCCATGCACGACAG  | 2B  | /      | /        | 1, 2 |

**Table 2:** Recommended primers for amplification and sequencing of the VP1 region of serotypes SAT1, 2 and 3.

## 5. SAMPLES PREPARATION

If the viral RNAs are not already extracted, extract total RNA according to your usual procedure. As an example, the following strains can be used as positive controls, depending on the targeted serotypes:

| Serotype | Positive control strain |
|----------|-------------------------|
| O        | O Manisa                |
| A        | A 22                    |
| C        | C Noville               |
| Asia 1   | Asia 1 Shamir           |
| SAT 1    | SAT 1 Botswana          |
| SAT 2    | SAT 2 Zimbabwe          |
| SAT 3    | SAT3/ZIM/1981           |

**Table 3:** Identification of positive controls to be used according to the targeted serotypes.

## 6. OPERATING PROCEDURE

### 8.1 Mix preparation

Prepare the mix for all reactions to be performed, including a NTC (Non Template Control), a positive control and for 1 or 2 additional reactions. Final volume is 25 µL per reaction, containing, for each selected primer pair, the following reagents:

| Reagents  | Final concentration | Volume (µL) |
|---|---------------------|-------------|
| DNase RNase free water                          | /                   | 12,6        |
| One-Step RT-PCR 5 X buffer                      | 1 X                 | 5,0         |
| dNTPs (10 mM each)                              | 400 µM              | 1,0         |
| Forward primer (100 µM)                         | 0,8 µM              | 0,2         |
| Reverse primer (100 µM)                         | 0,8 µM              | 0,2         |
| One-Step RT-PCR Enzyme (1 U / µL)               | 1 U                 | 1,0         |
| <b>Total to be dispensed per well (µL)</b>      | <b>20</b>           |             |
| <b>RNA volume to be dispensed per well (µL)</b> | <b>5</b>            |             |

**Table 4 :** Mix composition

Dispense 20 µL of mix into each tube and add 5 µL of RNA extract. Close the tubes and centrifuge briefly using a microcentrifuge.

## **8.2 RT-PCR thermal cycles**

The detailed programs, adapted to the different primer pairs, are described in Table 5. After assembly, put the tubes directly into the thermal cycler and start the program quickly.

|                                   | Type O               |             | Types A, C et Asia1  |             | Types SAT1-2-3       |             |
|-----------------------------------|----------------------|-------------|----------------------|-------------|----------------------|-------------|
| <i>Reverse Transcription (RT)</i> | 30 min at 50°C       |             | 30 min at 50°C       |             | 30 min at 50°C       |             |
| <i>Hot Start Taq activation</i>   | 10 min at 95°C       |             | 10 min at 95°C       |             | 10 min at 95°C       |             |
| <i>Denaturation</i>               | 1 min at 95°C        | x 35 cycles | 1 min à 95°C         | x 35 cycles | 1 min at 95°C        | x 35 cycles |
| <i>Hybridization</i>              | <b>1 min at 60°C</b> |             | <b>1 min at 55°C</b> |             | <b>1 min at 50°C</b> |             |
| <i>Elongation</i>                 | 1.5 min at 72°C      |             | 1.5 min at 72°C      |             | 1.5 min at 72°C      |             |
| <i>Final elongation</i>           | 5 min at 72°C        |             | 5 min at 72°C        |             | 5 min at 72°C        |             |

**Table 5:** Thermal cycles used for the amplification of the VP1 region

*Rq: It is possible to use other primer pairs referenced in paragraph 9 and/or in the literature or designed in-house. In this case, it is necessary to adapt the RT-PCR program applied to the hybridization temperatures of the primers, and the elongation time to the expected size of the amplified fragment.*

## **7. VISUALIZATION OF RESULTS**

Prepare a 2% agarose gel in 1 X TAE. For each sample, deposit 5 µL of amplified DNA mixed with 1 µL of loading buffer. For each migration, add at least one molecular weight marker. Allow 50 minutes of migration at 120 V (maximum amperage). Agarose gel migration conditions may vary depending on agarose concentration, gel size and amplified DNA size.

Immerse the gel in the ethidium bromide bath for 10-20 minutes. Visualize the amplified DNA under UV light using the imager. Save and export the gel pictures in jpeg format.

## **8. EXPRESSION OF RESULTS**

### **8.1. TEST VALIDATION**

The assay is validated if:

- No band is visible for negative RT-PCR controls (NTC)
- A band of the expected size is visible for positive RT-PCR controls

### **8.2. EXPRESSION OF RESULTS**

The expression of the results depends on the bands obtained:

- **Presence of a band at the expected size:** specific amplification. The DNA can be used for sequencing (in-house or outsourced).
- **No visible band:** negative result. Test other primers or other amplification conditions.
- **Presence of a band not corresponding to the expected size in addition to the specific band:** non-specific amplification. The DNA can be used for sequencing but it is preferable to test other primers or other amplification conditions.

## 9. ADDITIONAL PRIMERS

The primers described in Table 6 are called additional primers. They can be used to amplify VP1 or to sequence the amplification product obtained with the primers previously mentioned.

| Serotype | Sens | Name         | Sequence (5'-3')           | Region | Ref.    |
|----------|------|--------------|----------------------------|--------|---------|
| Pan FMD  | -    | NK72         | GAAGGGCCCCAGGGTTGGACTC     | 2A/2B  | 1, 2, 3 |
| Pan FMD  | -    | NK61         | GACATGTCTCTCTGCATCTG       | 2B     | 3       |
| O        | +    | O-CRH2F      | GAYTACGCSTACACSGCGTC       | VP3    | 1,2     |
| O        | +    | O-1D293F     | TGGAYAACACCACYAAAYCCAAC    | VP1    | 1,2     |
| O        | +    | O-1D296F     | ACAACACCACCAACCCAAC        | VP1    | 1,2     |
| O        | +    | O-1D296aF    | ATAACACCACTAATCCAAC        | VP1    | 2       |
| O        | +    | O-1D296bF    | ACAACACCACCAATCCAAC        | VP1    | 1,2     |
| O        | +    | O-1D296cF    | ATAACACCACCAATCCAAC        | VP1    | 2       |
| O        | +    | O-1C605aF    | TGGCTAGTGCTGGTAAAGACTTTGAG | VP3    | 2       |
| O        | +    | O-1C605bF    | TGGCTAGTGCCGGCAAGGACTTTGAG | VP3    | 2       |
| O        | +    | O-1C605cF    | TGGCTAGCGCCGGCAAGGACTTTGAG | VP3    | 2       |
| O        | +    | O-1C605dF    | TGGCTAGCGCCGAAAGGACTTTGAG  | VP3    | 2       |
| O        | -    | O-1D628bR    | GTTGGGTTGGTGGTGT           | VP1    | 2       |
| O        | -    | O-1D628R     | GTTGGGTTGGTGGTGGTTGT       | VP1    | 1,2     |
| O        | -    | O-1D628aR    | GTTGGATTAGTGGTGT           | VP1    | 1,2     |
| O        | +    | O-1C564F     | AATTACACATGGCAAGGCCGACGG   | VP3    | 6       |
| A        | -    | A-1D523R     | CGTTTCATRCGCACRAGRA        | VP1    | 1       |
| A        | +    | A-1D202aF    | TCAGCCACCTACTATTTCTCTGA    | VP1    | 2       |
| A        | +    | A-1D202bF    | GCAGCAACATACTACTTCTCTGA    | VP1    | 2       |
| A        | +    | A-1D202cF    | GCAGCAACCTACTATTTCTCTGA    | VP1    | 2       |
| A        | +    | A-1D202dF    | GCGGCCACTTACTACTTCTCTGA    | VP1    | 2       |
| A        | +    | A-1D202eF    | GCGGCCACCTACTATTTTCTGA     | VP1    | 2       |
| A        | +    | A-1D202fF    | GCGGCCACCTACTACTTTTCTGA    | VP1    | 2       |
| A        | +    | A-1D205F     | GCNACNTACTAYTTYTC          | VP1    | 5       |
| A        | -    | A-1D478aR    | CAGTGCTCCGTAGTTAAAGGATGA   | VP1    | 2       |
| A        | -    | A-1D478bR    | AATTGCACCGTAATTGAAGGATGC   | VP1    | 1,2     |
| A        | -    | A-1D478cR    | GAGTGCACCATAGTTGAAAGACGC   | VP1    | 2       |
| A        | -    | A-1D478dR    | AATCGCACCAAAGTTGAAGGAAGT   | VP1    | 2       |
| A        | -    | A-1D478eR    | AACTGCGCCGTAGTTGAAGGAGGC   | VP1    | 2       |
| A        | -    | A-1D478fR    | AATTGCGCCGTAGTTGAAGGATGC   | VP1    | 2       |
| C        | -    | C-1D535R     | ARAGYTCIGCICGYTTCAT        | VP1    | 1,2     |
| Asia     | +    | As1-1D205F   | GCRACGTACTACTTYTCRGACCT    | VP1    | 2,3     |
| Asia     | -    | As1-1D370R   | GTTGTAYACTGYGCCAGCACACG    | VP1    | 2,3     |
| SAT 1    | +    | SAT1-1D200F  | TGCGYGCGCCACGTACTAYTTCTC   | VP1    | 1,2     |
| SAT 1    | +    | SAT1-1D200aF | TGCGTGCGCCACGTATTATTCTC    | VP1    | 2       |
| SAT 1    | +    | SAT1-1D200bF | TGCGGCTGCTACGTACTACTTCTC   | VP1    | 2       |
| SAT 1    | +    | SAT1-1D200cF | TGCGYGCGCCACGTACTAYTTCTC   | VP1    | 2       |
| SAT 1    | +    | SAT1-1D200dF | TGCGTGCTTCCACGTACTACTTCTC  | VP1    | 2       |
| SAT 1    | +    | SAT1-1D200eF | TGCGYGCGCCACGTACTACTTCTC   | VP1    | 2       |
| SAT 1    | -    | SAT1-1D394R  | GGYTTGTA TTRCARTCACCGTTGTA | VP1    | 1,2     |

|       |   |               |                              |     |     |
|-------|---|---------------|------------------------------|-----|-----|
| SAT 1 | - | SAT1-1D394aR  | GGTTTGTAYTTGCAGTYGCCRTTGTA   | VP1 | 2   |
| SAT 1 | - | SAT1-1D394bR  | GGCTTGTACTTACAGTCACCATTGTA   | VP1 | 2   |
| SAT 1 | - | SAT1-1D394cR  | GGTTTGTAYTTGCAGTTGCCRTTGTA   | VP1 | 2   |
| SAT 1 | - | SAT1-1D394dR  | GGCYTGTACTTRCAGTCACCGTTGTA   | VP1 | 2   |
| SAT 1 | - | SAT1-1D394eR  | GGTTTGTAYTTGCARTCACCGTTGTA   | VP1 | 2   |
| SAT 2 | + | SAT2-P1-1223F | TGAACTACCACTTCATGTACACAG     | VP3 | 1,2 |
| SAT 2 | - | SAT2-D        | GGTGCGCCGTTGGGTTGCCA         | VP1 | 1,2 |
| SAT 2 | + | SAT2-1D209aF  | CCACTTACTACTTTTGTGACCTTGA    | VP1 | 2   |
| SAT 2 | + | SAT2-1D209bF  | CCACCTACTACTTTTGTGACCTTGA    | VP1 | 2   |
| SAT 2 | + | SAT2-1D209cF  | CCACCTACTATTTCTGTGACCTGGA    | VP1 | 1,2 |
| SAT 2 | + | SAT2-1D209dF  | CCACGTACTACTTCTGTGACCTGGA    | VP1 | 2   |
| SAT 2 | + | SAT-1D209F    | CCACATACTACTTTTGTGACCTGGA    | VP3 | 4   |
| SAT 2 | + | SAT2-1C513aF  | CACACACACAGACACACCCGCGATGGC  | VP3 | 2   |
| SAT 2 | + | SAT2-1C513bF  | CACACACACTGACACACCTGCGATGGC  | VP3 | 2   |
| SAT 2 | + | SAT2-1C513cF  | CACCCACACAGACACACCCGGCCATGGC | VP3 | 2   |
| SAT 2 | + | SAT2-1C513dF  | CACGCACACAGACACCCCGGCCATGGC  | VP3 | 2   |
| SAT 2 | + | SAT2-1C513dG  | CACGCACACGGACTCCCGCGATGGC    | VP3 | 2   |
| SAT 2 | + | SAT2-1C513dH  | CTCTCACACGGACTCGCGGAAGGC     | VP3 | 2   |

**Table 6:** Additional primers used for sequencing or amplification of the VP1 region of FMD virus.