



TECHNICAL EVALUATION OF TWO RAPID KITS FOR RABIES DIAGNOSIS

A. SERVAT & E. ROBARDET

EURL Workshop on Rabies, 15-16 June 2022

CONNAÎTRE, ÉVALUER, PROTÉGER

Objective

Assessment of the technical performance of two commercially available rapid immunochromatographic tests for rabies diagnosis, and to compare it to the Fluorescent Antibody Test

Characteristics :

- immunochromatographic tests = Lateral Flow Assays (LFAs) = Lateral Flow Device (LFDs)
- generally based on colloidal gold conjugated monoclonal antibodies that capture the antigen contained in a sample
- The antigen–antibody complex migrates on a nitrocellulose membrane and binds to a detection antibody fixed in the test zone “T” revealing a coloured line for a positive sample.
- Simple, rapid (results obtained in about 5 to 10 minutes) and that they do not request any special equipment

Samples

Sample status	Country of origin	Infection origin	Number of samples
Rabies positive	Morocco	natural	20
	Ukraine	natural	9
	Hungary	natural	4
	Serbia	natural	6
	Laboratory*	experimental	7
Rabies Negative	France	-	45

** Brains collected from mice experimentally infected with RABV or other lyssaviruses, and subsequently lyophilized*

91 brain samples were tested (46 rabies + and 45 rabies -)

LFD kits

- Petscreen *Canine Rabies antigen test*[®], from Global DX Ltd, United Kingdom,
- Redtest *rabies virus antigen rapid test*[®], Sigmed, Poland.



Intended use:

Petscreen = saliva and brain homogenates

Redtest = saliva and cerebrospinal fluid (=> brain specimens)

Results : FAT vs Global DX LFD

		FAT results	
		Positive	Negative
LFD results	Positive	3	0
	Negative	43	45

Specificity of **100%** (95% CI: 92.1 - 100%)

Sensitivity of only **6.5%** (95% CI: 1.37 - 17.9%)

Only 3 FAT positive samples (2 from Morocco and on freeze-dried sample infected with ABLV) were confirmed positive with this LFD

Results : FAT vs REDTEST LFD

		FAT results	
		Positive	Negative
LFD results	Positive	37	0
	Negative	9	45

Specificity of **100%** (95% CI: 92.1 - 100%)

Sensitivity of only **80,4%** (95% CI: 66.1 - 90.6%)

Redtest did not succeed to detect rabies antigens from the samples experimentally infected with Duvenhage, EBLV-1a, EBLV-1b, EBLV-2 and BBLV

Agreement between FAT and LFD results (qualitative results)

For each test, comparison between FAT and LFD (Kappa statistic test) => Determination of a Cohen's Kappa values

- ➔ FAT vs Global DX LFD results = poor agreement (Cohen's kappa = 0.065 - CI: 0 – 0.136) between the two tests.
- ➔ FAT vs Redtest LFD results = a good agreement (Cohen's kappa = 0.802 - CI: 0.682 – 0.923) between the two tests.

Kappa score	Interpretation
< 0	No agreement
0.0 — 0.20	Slight agreement
0.21 — 0.40	Fair agreement
0.41 — 0.60	Moderate agreement
0.61 — 0.80	Substantial agreement
0.81 — 1.00	Almost perfect agreement


LFD evaluations for rabies diagnosis

> 20 publications/reports from 2007 – 2021
> 1870 brain samples analyzed (LFD vs DFA), collections from about 20 countries
13 different LFD kits evaluated (*Anigen test from Bionote part of 17 studies*)

Performance :

Specificity generally equal to 100% (except 4 studies with sp between 93% and 99%).

Sensitivity quite heterogeneous according to kits



LFD for rabies diagnosis evaluations

8 kits with sensitivities ranging from 0 to 20%

Ubio - 0% (Eggerbauer et al. 2016)

Quicking - 0 to 6,6% (Eggerbauer et al. 2016)

Biogen - 0 to 2,2% (Eggerbauer et al. 2016)

Elabscience - 0 to 20% (Kimitsutmi et al. 2020, Klein et al. 2020)

Intermedical - 3% (Klein et al. 2020)

Lillitest - 1% (Klein et al. 2020)

Span Biotec - 0% (Klein et al. 2020)

Global DX - 6,5% (Servat & Robardet, unpub, 2020)

2 kits with sensitivities ranging from 22 to 70%

Creative Diagnostic - 22,7 to 70% (Eggerbauer et al. 2016)

GreenSpring - 31,8 to 50% (Eggerbauer et al. 2016)

13 different LFD evaluated

3 kits with sensitivities ranging from 80 to 100%

ADTEC - 88 to 94,3% (Kimitsutmi et al. 2020, Managitt et al. 2021)


Redtest - 80,4% (Servat & Robardet, unpub, 2020)

Anigen - 88,3 to 100% (except in 2 studies : 62% & 16% to 100%)



Conclusion

- ✓ Petscreen®: barely unable to detect rabies antigens in positive samples (93.5% of false negative) + poor agreement with FAT
- ✓ Redtest® : good performance (sp, se, agreement with FAT) on RABV infected brain specimens of various animal species, but failed to detect most of all other non-RABV lyssaviruses.
- ✓ Confirmation of the great heterogeneity of LFD performance.
 - Need for improvement of rapid tests
 - Need for batch to batch consistency
 - Need for more transparency of reagents used by manufacturers (Mabs etc.)
 - Need for qualification of kits before being used/authorized



Acknowledgements

This study was funded by the European Commission.

ANSES technical staff: V. Brogat, S. Kempff & E. Litaize

Laboratories that kindly supplied rabies positive brain materials :

- ✓ Drs Rudoï Oleksii & Ivan Polupan - *State Scientific and Research Institute of Laboratory Diagnostics and Veterinary and Sanitary Expertise, Kiev, Ukraine.*
- ✓ Dr Sami Darkaoui - *ONSSA, Rabat, Morocco.*
- ✓ Dr Melinda Kocsis and Dr Péter Malik - *National Food Chain Safety Office, Budapest, Hungary.*
- ✓ Drs Jelena Maksimovic Zoric & Vesna Milicevic - *Institute of Veterinary Medicine of Serbia, Belgrade, Serbia.*

Thank you for your attention !

