



**French Agency for Food,
Environmental &
Occupational Health
Safety**

Maisons-Alfort

**LABORATOIRE DE SANTE
ANIMALE**

ANIMAL HEALTH LABORATORY

Unité Zoonoses Bactériennes
Bacterial Zoonoses Unit

Laboratoire National de
Référence pour la
Brucellose, la Chlamydie,
la Fièvre Charbonneuse, la
Morve, la Tuberculose & la
Tularémie animales

*Animal Anthrax, Brucellosis,
Chlamydiosis, Glanders,
Tuberculosis & Tularemia
National Reference
Laboratory*

Laboratoire de Référence
O.I.E. / FAO
pour la Brucellose,
la Paratuberculose et
la Tuberculose bovine

*OIE /FAO
Reference Laboratory
for Brucellosis,
Bovine Paratuberculosis
and Tuberculosis*

Laboratoire de Référence de
l'UE pour la Brucellose

*EU Reference Laboratory for
Brucellosis*

Laboratoire de Référence de
l'UE pour les Maladies
Equines (Morve)

*EU Reference Laboratory for
Equine Diseases (Glanders)*

Supplier / Manufacturer	IDEXX Montpellier SAS
Commercial name	POURQUIER Brucellose Positive Control
Product reference	P00130
Batch number	428
Expiry date	08 October 2023
Description	Internal reference material (animal brucellosis serology) Vial of 1 mL.
Batch size	1437
Reference ANSES	D-18-5003/7928
Date of control	14 November 2018
Standard references	Annex C of 64/432/EC directive and OIE Manual of Diagnostic Tests and Vaccines for Terrestrial Animals (Chapt. 2.1.4) (<i>Current versions</i>)
Operating procedure	French standards AFNOR NF U 47-003, NF U 47-004, NF U 47-019 et NF 47-021

Control test	Methods	Results*
Rose Bengal Test	NF U 47-003	+++
Complement Fixation Test	NF U 47-004	186.67 IU/ml [93.33 ; 373.33]
Seroagglutination Test	NF 47-021	100 IU/ml [50 ; 200]
iELISA (serum individual)*	NF U 47-019	277 % [cut-off 110-120]

*The results are for indication only

Maisons-Alfort, 15 November 2018

Acacia FERREIRA VICENTE
UE / OIE / FAO Reference
Laboratory for Brucellosis

This report applies only to the items submitted to control.
Reproduction of this report is authorized in its integral form only.
This report includes 0 annex (es)

*kit iELISA Brucellosis serum Bovine batch 8020; short incubation; cut-off: 110-120%

